

Release Notes

Viedoc 4.15

Available for users on
2015-06-04

Release Notes

Viedoc 4.15

New functionality



List of new functionality that has been added to this release.

Link to eLearning

Admin and Designer now contains an eLearning link in the footer of each page. The link will open Viedoc User Guide in a new tab on your browser. The User Guide is continuously updated to reflect the new features of Viedoc 4, including Video tutorials and FAQs.

Export preview and graphs

The Export Preview has been reworked and enhanced with new graph views. The data selection of Export remain the same but when clicking Export Preview the table and graphs that are shown are new. There are three different graphs available; pie, column, and line.

Alerts

It is now possible to define alerts within Viedoc Designer. An alert is based on a condition that when evaluated to true will trigger a notification to the user(s) having a specific role. The notification is based on the internal Viedoc 4 messaging system. Alerts will appear as messages in the Message view available on the Landing page. It is also possible to choose the email option in the alert setup, this will send the alert notification as an email.

OID and Labels

An overview page with editing possibilities has been added to handle output IDs and output labels. The page will contain all items ordered by item groups and forms. Each item is displayed with its Field label, Field ID, Output field label and Output field ID. It is also possible to search for specific forms in order to find specific items in the overview page.

Release Notes

Viedoc 4.15

Updated functionality and bug fixes



List of functionality that has been updated in this release.

Updated Functionality

- Study settings indicators.
- Export includes all forms, with or without data.
- Configuration report includes Validation result, SDV, and Alerts.
- Validation rules updated, now requires at least one active Role.
- EventDate edit checks can now be made cross form
- Possibility to edit and delete a role from inside the Role edit page is removed, the functions are still available at Roles overview page.
- Export settings have “Group data by form” selected by default, this improves the loading of preview significantly.

Bug Fixes

The following corrections have been implemented in this release of Viedoc:

- Avatar settings can now be made, changed, and cancelled in Selection View Settings.
- Viedoc Message view now works with IE8.
- Empty audit trail records are no longer created when applying revision.
- Fields with decimal values are now set when applying revision.
- Validation errors will only highlight the item group, not the radio and checkboxes inside.
- Number of validation errors now show correct when validating combined expressions.
- Randomization validation errors is now displayed correctly.
- Randomization definition is now included in validation.
- Expressions with _ (underscore) will not result in a validation error.
- Deleted items will not be included in validation.
- SAS export content will not include line breaks.
- System variables are no longer truncated in SAS.

Release Notes

Viedoc 4.15



List of updated functionality continues.

- Language import with form ID containing _ (underscore) is now possible.
- It is now possible to apply revision to a planned visit without study design.
- Audit trail records now show correct when a function returns null.
- Audit trail records no longer rounds the value when a function returns decimals.
- Patient card variables in the Selection view now show correct information for new and existing patients.
- Import of clinical data now includes all 5 variables of Selection View settings.
- Hidden fields are no longer shown in export preview.

Release Notes

Viedoc 4.15

Known limitations



List of limitations that have been found with this and earlier releases of Viedoc, and will be changed in future releases.

- Only one (selected) language can be imported from each CDISC ODM file
- Only one (selected) design can be imported from each CDISC ODM file
- When the study workflow is saved, the language in the design settings are lost and have to be reset
- Visibility conditions (entry conditions) cannot be set for activities in common events.
- The next and previous links may not correctly reflect the number of subject filtered in some instances. Going back to subject selection and re-selecting a patient solves the problem.
- Opening a form might fail in some instances when using IE 8.
- The PDFs generated as Empty CRFs will not display all code list items for radio and checkboxes if these have been configured with many code list items in a vertical layout.
- The code list tables for dropdowns are not displayed completely in the CRF PDF if the code list has many items.
- A field that is required for SDV but is hidden on a form (due to for example visibility conditions) cannot normally be marked as having been SDV:ed, but if the entire form is marked as having been SDV:ed then even hidden fields are included and marked as having been SDV:ed.
- When using the new patient table view some options default back to the patient card view, for example choosing the order of display, the site, or using the search capability. The user has to enter the table view again, but the option chosen is preserved and is then displayed correctly in the table view.
- The affected forms count summary shown in admin while applying a revision is implemented only for production site. Demo/Training site forms are not included in this summary count.
- Deletion of a common event is shown with the form name strike through instead of the summary format details.
- Selection View settings (former Patient Card settings) is limited to 5 variables. Old designs may have settings containing more than 5 variables. When importing such a design the variables 1-5 will be used and all other variable will be ignored.
- Validation of selection view settings and patient id generation settings for deleted items are not performed.
- Export preview with graph pages result in a Java script error on IE8.

Release Notes

Viedoc 4.15

Signoff



Signature Confirmation by
the product owner.

Uppsala, Sweden 2015-05-28

Henrik Blombergsson

PCG Solutions