Available for users on 2014-02-19

Updated functionality and bug fixes

List of new functionality that has been updated in this release.

- Added the possibility to import study clinical data from a CDISC ODM file (XML)
- Added the possibility to set and show one form item in the unscheduled events dashboard
- Added support for "SiteCode" to be part of the auto-generated subject identifier
- Improved support for Internet Explorer 8

The following limitation has been corrected in this release of Viedoc, and is now fixed:

• The item order in exports does not match the item order in forms

Known limitations

List of limitations that have been found with this and earlier releases of Viedoc, and will be changed in future releases.

Ц

- Queries cannot be exported
- The export preview feature is heavy to load and might cause the application to feel unresponsive
- Only one (selected) language can be imported from each CDISC ODM file
- Only one (selected) design can be imported from each CDISC ODM file
- When a study is exported in CDISC ODM format a study id is generated and added to the file exported. If this file is then imported into another study then the study id is not reset, and appears to be the same for both studies. This study id is not used internally in Viedoc so this has no effect on Viedoc, but can be confusing if export files are reviewed manually.

Signoff

Signature Confirmation by the product owner.

Uppsala, Sweden 2014-02-12

burgston en

Henrik Blombergsson

PCG Solutions